


CODE
INCUS 07.17

Audio-HA® Incus

HA: Bioactive calcium hydroxylapatite ceramic

Prosthesis for insertion between mobile malleus and stapes.
Its compact design allows a good view of the stapes head and of the handle of the malleus while positioning the prosthesis.
After measuring (Microsurgical Measuring Device: Catalog Number 02.14), the prosthesis can be trimmed using a diamond bur under continuous irrigation.



Dimensions

Total length (mm):


3.25 - 3.75

Columella outside diameter: 2.00 mm

Columella inside diameter: 1.20mm

Packaging

Dispenser box containing 6 sterile,
individually packaged implants.

	Total length (L) mm	Code
	3.25	INCUS 07.17.325
	3.75	INCUS 07.17.375

