

CODE  
SPL 03.20

## Titanium Piston

Titanium GR. 2

### Dimensions

Total length:


4.00 - 4.25 - 4.50 - 4.75

5.00 - 5.25 - 5.50 - 5.75

Piston diameter: 0.6 mm.

### Packaging

Dispenser box containing 6 sterile, individually packaged implants.

	Total length (L) mm	Code
	4.00	SPL 03.20.400
4.25	SPL 03.20.425	
4.50	SPL 03.20.450	
4.75	SPL 03.20.475	
5.00	SPL 03.20.500	
5.25	SPL 03.20.525	
5.50	SPL 03.20.550	
5.75	SPL 03.20.575	

CODE  
SPL 03.21

## Titanium Piston

Titanium GR. 2

### Dimensions

Total length:


4.00 - 4.25 - 4.50 - 4.75

5.00 - 5.25 - 5.50 - 5.75

Piston diameter: 0.4 mm.

### Packaging

Dispenser box containing 6 sterile, individually packaged implants.

	Total length (L) mm	Code
	4.00	SPL 03.21.400
4.25	SPL 03.21.425	
4.50	SPL 03.21.450	
4.75	SPL 03.21.475	
5.00	SPL 03.21.500	
5.25	SPL 03.21.525	
5.50	SPL 03.21.550	
5.75	SPL 03.21.575	

CODE  
SPL 03.22

## Titanium Piston


Titanium GR. 2

**Dimensions**

Total length:  
4.00 - 4.25 - 4.50 - 4.75  
5.00 - 5.25 - 5.50 - 5.75  
Piston diameter: 0.5 mm.

**Packaging**

Dispenser box containing 6 sterile, individually packaged implants.

	Total length (L) mm	Code
	4.00	SPL 03.22.400
	4.25	SPL 03.22.425
	4.50	SPL 03.22.450
	4.75	SPL 03.22.475
	5.00	SPL 03.22.500
	5.25	SPL 03.22.525
	5.50	SPL 03.22.550
	5.75	SPL 03.22.575

CODE  
SPL 03.23

## Titanium Piston


Titanium GR. 2

**Dimensions**

Total length:  
4.00 - 4.25 - 4.50 - 4.75  
5.00 - 5.25 - 5.50 - 5.75  
Piston diameter: 0.8 mm.

**Packaging**

Dispenser box containing 6 sterile, individually packaged implants.

	Total length (L) mm	Code
	4.00	SPL 03.23.400
	4.25	SPL 03.23.425
	4.50	SPL 03.23.450
	4.75	SPL 03.23.475
	5.00	SPL 03.23.500
	5.25	SPL 03.23.525
	5.50	SPL 03.23.550
	5.75	SPL 03.23.575